



ST. LUKE'S

Glucose Meter Test

1. The patient test result displays as ▲ 423mg/dL. What does this mean?
 - a. It is a high critical result
 - b. It is greater than 423
 - c. It is an invalid result
2. Which patients are considered critically ill by SLH definition?
 - a. CVU, post CABG patient, Level 1
 - b. ICC Back, DKA insulin drip and arterial line, Level 1
 - c. Overdose patient on a vent, Level 1
 - d. ICC Front, trach/vent patient with hospice care and no lines, Level 2
3. When you open a new box of quality controls you must label each bottle with the expiration date. They expire ___ days from the date you open them
 - a. 30
 - b. 60
 - c. 90
4. Why should the 1st drop of blood be wiped from the patient's finger with gauze and the test be performed on the 2nd drop?
 - a. To assure the finger is cleaned properly
 - b. So tissue fluid does not alter the blood sample
 - c. Because the first drop of blood is too thick for the sample strip
5. Accurate bedside glucose results rely on good circulation to the fingers. If a patient has poor circulation in the fingers or is admitted as a Level 1 patient, DO NOT perform a bedside glucose test.
 - a. True
 - b. False
6. Within a few seconds after applying a drop of blood to the test strip, the FSPP meter should beep and indicate "SAMPLE ACCEPTED." If this message does not appear, how long do you have to apply a second drop of blood?
 - a. 20 seconds
 - b. 30 seconds
 - c. A second drop of blood cannot be applied
7. Why is it important to keep the glucose meter flat (horizontal) when running controls or a patient test?
 - a. You will get an error code if you do not keep the glucose meter flat
 - b. Holding the glucose meter upright causes liquids to enter the test port which will require the meter to be replaced
 - c. The glucose meter will turn off if tipped up or down



ST. LUKE'S

Glucose Meter Test

8. What steps would you take if the Quality Control FAILS?
 - a. Check the quality control bottles and make sure they have not expired. Replace if expired
 - b. Check to make sure the correct control (LO or HI) was analyzed
 - c. Remix the control, expel any air bubbles, and re-run the control
 - d. All of the above
9. Comment Code "4" is entered with the patient test result to:
 - a. Acknowledge the result is a critical value
 - b. Delete a patient test when you make an error
 - c. Indicate you repeated the test
10. If you made an error running a patient glucose test and did **not** enter Comment Code #4 when the result was displayed, how would you correct the record?
 - a. Call the laboratory and report your error
 - b. Complete a POC Corrected Report Form and forward as instructed
 - c. Call the IT Service Desk for help
11. If the bar-code scanner fails to scan, what would you do to troubleshoot? (select all that apply)
 - a. Replace the two AA batteries if battery indicator is low
 - b. Tap the glucose meter on the countertop to "reset" scanner and batteries
 - c. Clean the laser scanning window with a soft cloth
12. Your patient test result displays as >500mg/dL. What does this mean? (select all that apply)
 - a. The result is too high for the FSPP Pro meter to read
 - b. The result is less than 500mg/dL
 - c. A laboratory glucose test should be obtained to get an accurate glucose result.
13. What number do you enter into the FSPP meter for the Operator ID?
 - a. Your employee ID number, which can be scanned into the meter using the barcode on your ID badge
 - b. Your Social Security number
 - c. All 9's (999999999)
14. What patient number do you enter into the FSPP meter for the PATIENT ID?
 - a. Social Security number
 - b. Current FIN
 - c. Medical Record Number (MR)
15. How often must you run both levels of quality controls?
 - a. Once a month
 - b. Once a week
 - c. Every day of patient testing



ST. LUKE'S

Glucose Meter Test

16. If your FSPP meter is malfunctioning, where do you take it to obtain a replacement?
 - a. TriMedx
 - b. IT
 - c. The Laboratory (24 hours a day, 7 days a week)
17. If a patient has not been assigned a FIN and a glucose is needed emergently (newborns, ER urgent situation, etc.) which steps would you do?
 - a. Perform the testing using all 9's as the patient ID
 - b. Complete the POC Corrected Report Form including the patient information once registration is complete
 - c. Send completed POC Corrected Report Form to the lab for follow up
 - d. All of the above
18. Why is it important to make sure the patient ID band is replaced when a patient is transferred in from another facility?
 - a. The barcode on the ID band from another facility can be scanned, but will be the wrong FIN and test results will not be transmitted
 - b. As long as the barcode on the ID band scans, it is acceptable to use for testing
 - c. The barcode from another hospital will not scan
19. As a matter of good clinical practice, caution is advised in the interpretation of NEONATAL glucose values below _____. Order a lab glucose for more accurate results.
 - a. 50 mg/dL
 - b. 100 mg/dL
 - c. 10 mg/dL
20. Which type of patient(s) may be used for glucose testing on the FSPP glucose meter:
 - a. FBC, Birthing Mom
 - b. ER, Teenager with a broken arm
 - c. 2E, Elderly patient with pneumonia
 - d. 3W, Middle aged man with a total knee replacement
 - e. All of the above
21. You are caring for a Level 1 critically ill patient. Which of the following is the correct method for obtaining a blood glucose result (choose all that apply)?
 - a. Obtain a point of care capillary glucose sample
 - b. Order a glucose level
 - c. Obtain venous blood sample and send to lab in a green tube
 - d. All of the above