



CERNER FLASH

6/7/2018

BRIDGE BLOOD TRANSFUSION NOW WORKING

Tag issue resolved. Please use Bridge for all Blood Transfusions going forward. Thanks!

Standard Transfusion Steps:

1. Access Bridge through the button on the toolbar within Powerchart/Surginet
2. Scan patient's wristband to access patient's chart in Bridge
3. Click 'Start Transfusion'
4. Follow the Prechecks
5. Scan the recipient tag on the blood
6. Scan the 4 barcodes on the blood bag in the following order:
 - a. Unit number
 - b. Blood product
 - c. Expiration date
 - d. Donor blood type
7. Select START after spiking the blood and starting the transfusion