

FetaLink Launch from mPages Issue Resolved

Prior Issue:

☑ FetaLink was intended to launch from PowerChart to allow for review of strips that directly correlated. However “Launch” button did not function.

Resolution Effective Immediately:

☑ FetaLink launches from PowerChart as expected.

1. From desired mPage, click **Fetal Monitoring**.
2. If patient has been seen in any McLaren office or hospital (with Cerner) their past fetal monitoring will be displayed based on date of monitoring. Note: Nothing will display if monitors are not appropriately associated/disassociated to patient or if patient has not been seen.
2. Check mark those that you wish to review; as you select the date; a Launch button will appear in the bottom right corner.
3. Once Launch is selected, FetaLink** will open to designated date of service.

****Note: This will ONLY work on computers with FetaLink installed on them.**

The screenshot shows the PowerChart interface for a patient's Women's Health View. The main content area displays a table of fetal monitoring records. The table has the following columns: Start Date/Time, Stop Date/Time, EGA at Start, and Reason for Monitoring. A single record is shown for June 16, 2019, 19:19, with a stop date of June 17, 2019, 05:14, and a reason for monitoring of labor. A checkmark is visible in the left margin of the record, and a yellow 'Launch' button is highlighted in the bottom right corner of the table.

Start Date/Time	Stop Date/Time	EGA at Start	Reason for Monitoring
June 16, 2019 19:19	June 17, 2019 05:14	40 Weeks 4 Days	labor

